

## Recombinant Spike protein S[679-833]

PX-COV-P011

### OVERVIEW


---

<b>ORIGIN SPECIES</b>	SARS-CoV2
<b>FRAGMENT</b>	Partial [679-833]
<b>PROTEIN DELIVERED WITH TAG</b>	N-terminal Histidine tag
<b>MOLECULAR WEIGHT WITH TAG IF ANY</b>	19.2kDa
<b>DELIVERY CONDITION</b>	Dry Ice

### PRODUCT INFORMATION

---

<b>EXPRESSION SYSTEM</b>	Prokaryotic expression
<b>HOST</b>	E.coli
<b>PURITY</b>	90%
<b>PROTEIN ACCESSION</b>	MN908947
<b>FORM</b>	Frozen
<b>BUFFER</b>	PBS, pH7.5, 4M urea
<b>STABILITY &amp; STORAGE</b>	ProteoGenix did not perform stability or storage tests for your specific end products. Standard storage advice: Store at 4°C for short term; Freeze and store at -20°C or -80°C for long term. Avoid freeze/thaw cycles



### More info

#### SYNONYMS

Spike protein fragment  
679-833