

Human HSP70 Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P1057
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Full-length
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	73,01 kDa
FORM	liquid
BUFFER	PBS, imidazole 300mM, pH7.5 without cleavage
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P0DMV8
ALIAS/SYNONYMS	HSP70, HSP70-1

