

Human UHRF1(400-660) Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P1062
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Partial
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	55,97 kDa
FORM	liquid
BUFFER	TrisHC 50mM, 10mM reduced glutathion pH8
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q96T88
ALIAS/SYNONYMS	ICBP[400-660], transcription factor ICBP90

