

Leishmania CPACter Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P1064
ORIGIN SPECIES	Leishmania
FRAGMENT	Partial
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	27,85 kDa
FORM	liquid
BUFFER	PBS, imidazole 300mM
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q9BIE1
ALIAS/SYNONYMS	CPACter (E. coli optimized), Cysteine Peptidase A C-terminal region, Cysteine proteinase-like protein

