

Human DJ1 Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P1084
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Full-length
TAG	No
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	47,02 kDa
FORM	liquid
BUFFER	TrisHC 50mM, NaCl 150mM, EDTA 1mM, DTT 1mM
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q99497
ALIAS/SYNONYMS	DJ1, protein DJ-1 (Alias: Oncogene DJ1 Parkinson disease protein 7)

