

Mouse ENO1 Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P1088
ORIGIN SPECIES	House Mouse
FRAGMENT	Full-length
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	50,68 kDa
FORM	Frozen
BUFFER	PBS, pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P17182
ALIAS/SYNONYMS	ENO1 Protein, Alpha Enolase Protein, Alpha enolase like 1 Protein, ENO1L1 Protein, Enolase 1 Protein, MBP-1 Protein, MBPB1 Protein, MPB1 Protein, MYC promoter-binding protein 1 Protein, NNE Protein, Phosphopyruvate hydratase Protein, Alpha-enolase Protein

