

Human FGF2b Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P1089
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Partial
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	17,76 kDa
FORM	liquid
BUFFER	PBS, Urea 8M
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P09038
ALIAS/SYNONYMS	FGF2b, basic fibroblast Growth Factor proteins, FGF2, Fibroblast Growth Factor proteins 2, FGF-2, bFGF, Heparin-binding Growth Factor proteins 2, HBGF-2

