

Human GC_006 Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P1092
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Partial
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	39,70 kDa
FORM	liquid
BUFFER	PBS, imidazole 300mM
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q5C9Z4
ALIAS/SYNONYMS	GC_006 (NOM1[230-570]), nucleolar protein with MIF4G domain 1

