

Human HA-HSP27 Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P1109
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Full-length
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	25,86 kDa
FORM	liquid
BUFFER	PBS, imidazole 300mM. PBS only if dialysis
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P04792
ALIAS/SYNONYMS	HA-HSP27, Heat Shock 27kDa Protein 1 / heat shock protein beta-1

