

Human HSP110wt Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P1113
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Full-length
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	98,13 kDa
FORM	liquid
BUFFER	PBS, Urea 6M
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q92598
ALIAS/SYNONYMS	HSP110wt, heat shock protein 105 kDa isoform 1

