

Human UHRF1-ICBP90 tagged Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P1115
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Full-length
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	91,09 kDa
FORM	liquid
BUFFER	PBS, imidazole 300mM
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q96T88
ALIAS/SYNONYMS	ICBP90, transcription factor ICBP90, E3 ubiquitin-protein ligase UHRF1, Inverted CCAAT box-binding protein of 90 kDa, Nuclear protein 95, Nuclear zinc finger protein Np95, HuNp95, hNp95, RING finger protein 106, RING-type E3 ubiquitin transferase UHRF1, T

