

Human KELCH Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P1118
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Partial
TAG	His
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	29 kDa
FORM	liquid
BUFFER	PBS, Urea 4M
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q8WZ60
ALIAS/SYNONYMS	KELCH, kelch-like 6

