

Human MSH5Cter Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P1126
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Partial
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	49,54 kDa
FORM	liquid
BUFFER	NaH ₂ PO ₄ 100mM, TrisHCl 10mM, Urea 8M
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	O43196
ALIAS/SYNONYMS	MSH5Cter, mutS protein homolog 5 isoform c

