

Human NDKA Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P1130
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Full-length
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	18,82 kDa
FORM	liquid
BUFFER	PBS, imidazole 250mM, pH 8.0
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P15531
ALIAS/SYNONYMS	NDKA, nucleoside diphosphate kinase A

