

Human NPC1L1 Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P1131
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Partial
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	20,81 kDa
FORM	liquid
BUFFER	PBS, imidazole 300mM in native conditions.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q9UHC9
ALIAS/SYNONYMS	NPC1L1 [450-620], NPC1L1 protein, truncated Niemann-Pick C1-like protein 1

