

Xenopus rtTA Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P1148
ORIGIN SPECIES	Xenopus
FRAGMENT	Full-length
TAG	
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	30,56 kDa
FORM	liquid
BUFFER	PBS, imidazole 400mM
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	ABC65842.1
ALIAS/SYNONYMS	rtTA, rtTA-M2.alt [Expression vector peTetOn(PacI)], doxycycline-activated transcriptional activator

