

Human TIVAMP Human Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P1161
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Partial
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	26,11 kDa
FORM	liquid
BUFFER	NaH ₂ PO ₄ 100mM, TrisHCl 10mM, Urea 8M
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P51809
ALIAS/SYNONYMS	TIVAMP Human, vesicle-associated membrane protein 7 isoform 1

