

Human TRAP1 Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P1164
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Full-length
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	82,16 kDa
FORM	liquid
BUFFER	PBS, Urea 8M, imidazole 400mM
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q12931
ALIAS/SYNONYMS	TRAP1, TNF Receptor-associated Protein 1 (TRAP1), heat shock protein 75 kDa

