

Human UFD Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P1166
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Partial
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	30,24 kDa
FORM	liquid
BUFFER	75mM NaH ₂ PO ₄ , NaCl 150mM, Urea 4M
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q92890
ALIAS/SYNONYMS	UFD Δ Nter, ubiquitin fusion-degradation 1 like protein

