

IGFBP1 protein - Human IGFBP1 full length Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P2077
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Full-length
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	28.74kDa
FORM	Lyophilized
BUFFER	PBS, pH 7.5 with 0.2mM β -Met and glycerol 2.5%
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P08833
ALIAS/SYNONYMS	IGFBP1 full length, insulin-like Growth Factor proteins binding protein 1, Insulin-like Growth Factor proteins-binding protein 1