

Human IGFBP1 Nter Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P2078
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Partial
TAG	
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	17.59 kD
FORM	Lyophilized
BUFFER	PBS, pH 7.5 with 0.2mM β -Met
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P08833
ALIAS/SYNONYMS	IGFBP1 Nter, insulin-like Growth Factor proteins binding protein 1, Insulin-like Growth Factor proteins-binding protein 1

