

Human PAD4 arginine deiminase type-4 Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P2084
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Pro663
TAG	N-terminal GstTag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Insect
MOLECULAR WEIGHT	100.04 kDa
FORM	
BUFFER	PBS, pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q9UM07
ALIAS/SYNONYMS	Protein-arginine deiminase type-4 Human, peptidyl arginine deiminase, type IV, Protein-arginine deiminase type-4 ou HL-60 PAD ou Peptidylarginine deiminase IV ou Protein-arginine deiminase type IV, PAD4