

VSVG Protein- Vesicular stomatitis virus VSVG IND1 Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P2100
ORIGIN SPECIES	Vesicular stomatitis virus
FRAGMENT	Full-length
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Insect
MOLECULAR WEIGHT	58.44 kDa
FORM	Frozen
BUFFER	Tris-HC 50mM, NaCl 150mM, pH7.5
STABILITY & STORAGE	At 4°C for short term period (less than a week) At -20°C or -80°C for long term period It is important to avoid avoid freezing/thawing cycles 20-40% glycerol may be added to improve cryoprotection
UNIPROT	Q89015
ALIAS/SYNONYMS	VSVG IND1, G protein ou Glycoprotein