

Human Serum amyloid A-1 protein(SAA1) Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P3003
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Partial
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	13.2 kDa
FORM	Frozen
BUFFER	Tris-HC 50mM, NaCl 150mM, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P0DJ18
ALIAS/SYNONYMS	Serum amyloid A-1 protein, SAA, Amyloid fibril protein AA