

Rabies virus Nucleoprotein Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P3020
ORIGIN SPECIES	Rabies virus
FRAGMENT	Full-length
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Insect
MOLECULAR WEIGHT	51.58 kDa
FORM	Frozen
BUFFER	Tris-HC 50mM, NaCl 150mM, pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q0PNB8
ALIAS/SYNONYMS	Nucleoprotein

