

MousePAD2 Protein-arginine deiminase type-2 Mouse Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P3022
ORIGIN SPECIES	House Mouse
FRAGMENT	Full-length
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Insect
MOLECULAR WEIGHT	79.36 kDa
FORM	Lyophilized
BUFFER	50mM Tris-Cl pH 8.0, 150mM NaCl
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q08642
ALIAS/SYNONYMS	PAD-H19, PAD2, PDI2, Peptidylarginine deiminase II