



Human EGF Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P3041
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Partial
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	33.19 kDa
FORM	Frozen
BUFFER	PBS,pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P01133
ALIAS/SYNONYMS	Epidermal Growth Factor proteins, EGF, URG, Beta-Urogastrone, Pro-epidermal Growth Factor proteins

