

Human ICAM1, CD54 Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P3044
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Partial
TAG	
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	50.34 kDa
FORM	Frozen
BUFFER	PBS,pH 7.5 urea +8M
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P05362
ALIAS/SYNONYMS	CD54, BB2, P3.58, Human Rhinovirus Receptor, Major group rhinovirus receptor, Intercellular Adhesion Molecule 1, Major group rhinovirus receptor, ICAM-1