

Human PD CD6IP Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P3047
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Partial
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	44.48 kDa
FORM	Frozen
BUFFER	PBS,pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q8WUM4
ALIAS/SYNONYMS	PD CD6-IP, AIP1, Alix, DRIP4, HP95, ALG-2 Interacting Protein X,Programmed Cell Death Protein 6 Interacting Protein, PD CD6IP