

Human SOX2 Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P3053
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Unknown
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	22.96 kDa
FORM	Frozen
BUFFER	PBS,pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P48431
ALIAS/SYNONYMS	ANOP3, MCOPS3, SRY-Box 2, Transcription factor SOX-2, Sex Determining Region Y Box Protein 2, SOX2