

Human Protein S100 beta (S100B) recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P4007
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Full-length
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	11.78kDa
FORM	Lyophilized
BUFFER	PBS, pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P04271
ALIAS/SYNONYMS	S100-B, NEF, S100Beta, S-100 protein subunit beta