

Human Intra Acrosomal,SP-10 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P4029
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Full-length
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	29.7kDa
FORM	Lyophilized
BUFFER	PBS, pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P26436
ALIAS/SYNONYMS	Acrosomal protein SP 10, Acrosomal protein SP-10, Acrosomal vesicle protein 1, ACRV1, ASPX_HUMAN, D11S4365, Msa63, SP 10, Sp10, SPACA2, Sperm protein 10