

Human Golgi membrane protein 1 (GP73) recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P4051
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Partial
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	42.35kDa
FORM	Lyophilized
BUFFER	50 mM Tris-HCl pH 8, 150 mM NaCl
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q8NBJ4
ALIAS/SYNONYMS	GOLM1, GOLPH2, C9orf155, Golgi Membrane Protein 1, Golgi Phosphoprotein 2, GP73