

CD81 Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P4091
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Phe113~Lys201
TAG	N-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	10.59kDa
FORM	
BUFFER	PBS pH7.5
STABILITY & STORAGE	Store at - 20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage.
UNIPROT	P60033
ALIAS/SYNONYMS	TAPA1, TSPAN28