

CD317 Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P4118
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Ser50~Ser161
TAG	N-terminal Gst Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	39.49kDa
FORM	
BUFFER	PBS pH7.5
STABILITY & STORAGE	Store at - 20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage.
UNIPROT	Q10589
ALIAS/SYNONYMS	/

