

CD160 antigen(1-194)

OVERVIEW

PRODUCT REFERENCE	PX-P4287
ORIGIN SPECIES	
FRAGMENT	Met1-Ala194
TAG	
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	22.52kDa
FORM	
BUFFER	PBS, pH=7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	O95971
ALIAS/SYNONYMS	Natural killer cell receptor BY55,CD_antigen: CD160,BY55,CD160-TM

