

N-acetylneuraminate cytidyltransferase(neuA) His Tag

OVERVIEW

PRODUCT REFERENCE	PX-P4334
ORIGIN SPECIES	
FRAGMENT	Met1-Ser228
TAG	
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	26.22kDa
FORM	
BUFFER	PBS pH 7.4
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P0A0Z8
ALIAS/SYNONYMS	siaB, synB, CMP-N-acetylneuraminic acid synthase, CMP-NeuNAc synthase, CMP-sialic acid synthase

