

Mite allergen Blo t 5(BLOT5)

OVERVIEW

PRODUCT REFERENCE	PX-P4389
ORIGIN SPECIES	
FRAGMENT	Gln18-Gln134
TAG	
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	16.22KDa
FORM	
BUFFER	50 mM Tris-HCl pH 8.0, 150 mMNaCl
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	O96870
ALIAS/SYNONYMS	BLOT5/Mite allergen Blo t

