

Mucin-4(MUC4) EGF Like 1 and 2

OVERVIEW

PRODUCT REFERENCE	PX-P4396
ORIGIN SPECIES	
FRAGMENT	Gln1875-Glu2117
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	30.26KDa
FORM	
BUFFER	50 mM Tris-HCl pH 7.5, 150 mM NaCl, 5 mM GSH, 0.5 mM GSSG
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q99102
ALIAS/SYNONYMS	MUC4/Mucin-4, ASGP-1

