

Immunoglobulin G-binding protein G(spg)

OVERVIEW

PRODUCT REFERENCE	PX-P4551
ORIGIN SPECIES	
FRAGMENT	Lys314~Glu497
TAG	
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	20.96KDa (25-35kDa)
FORM	
BUFFER	PBS pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P19909
ALIAS/SYNONYMS	IgG-binding protein G

