

IL2 Protein - Interleukin-2 Protein (IL2)

OVERVIEW

PRODUCT REFERENCE	PX-P4557
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Ala21-Thr153
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	12-14kDa
FORM	
BUFFER	PBS pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P60568
ALIAS/SYNONYMS	T-cell Growth Factor proteins、TCGF、INN: Aldesleukin

