

Interleukin-4

OVERVIEW

PRODUCT REFERENCE	PX-P4558
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	His24-Ser153
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	16.20kDa
FORM	
BUFFER	PBS pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P05112
ALIAS/SYNONYMS	B-cell stimulatory factor 1、BSF-1、Binetrakin、Lymphocyte stimulatory factor 1、Pitrakinra

