

Interleukin-18(IL18)

OVERVIEW

PRODUCT REFERENCE	PX-P4562
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Tyr37-Asp193
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	19.45kDa
FORM	
BUFFER	PBS pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q14116
ALIAS/SYNONYMS	IGIF, IL1F4,Iloctadekin,Interferon gamma-inducing factor2 Publications,IFN-gamma-inducing factor2 Publications,Interleukin-1 gamma,IL-1 gamma

