

EGFP Protein - Enhanced green fluorescent protein(egfp)

OVERVIEW

PRODUCT REFERENCE	PX-P4577
ORIGIN SPECIES	Human cytomegalovirus (HHV-5) (Human herpesvirus 5)
FRAGMENT	Met1-Lys239
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	32kDa
FORM	
BUFFER	PBS pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	C5MKY7
ALIAS/SYNONYMS	