

Interleukin-1 alpha(IL1A)

OVERVIEW

PRODUCT REFERENCE	PX-P4588
ORIGIN SPECIES	
FRAGMENT	Ser113-Ala271
TAG	
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	19.2kDa
FORM	
BUFFER	50 mM Tris-HCl pH 8.0, 150 mM NaCl
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P01583
ALIAS/SYNONYMS	IL1F1,Hematopoietin-1

