

Keratin, type I cytoskeletal 19(KRT19)

OVERVIEW

PRODUCT REFERENCE	PX-P4594
ORIGIN SPECIES	
FRAGMENT	Met1-Leu400
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	45.1kDa
FORM	
BUFFER	0.02% Sarcosyl
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P08727
ALIAS/SYNONYMS	Cytokeratin-19, K19, KRT19,CK-19

