

Microtubule-associated protein tau(MAPT)(1-441)

OVERVIEW

PRODUCT REFERENCE	PX-P4601
ORIGIN SPECIES	
FRAGMENT	Met1-Leu441
TAG	
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	
FORM	
BUFFER	PBS pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P10636-8
ALIAS/SYNONYMS	Neurofibrillary tangle protein, Paired helical filament-tau, MAPTL, MTBT1, TAU

