

MCherry Protein - MCherry fluorescent protein(mCherry)

OVERVIEW

PRODUCT REFERENCE	PX-P4612
ORIGIN SPECIES	
FRAGMENT	Met1-Lys236
TAG	
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	
FORM	
BUFFER	PBS pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	X5DSL3
ALIAS/SYNONYMS	

