

Integrator complex subunit 3(Ints3)

OVERVIEW

PRODUCT REFERENCE	PX-P4683
ORIGIN SPECIES	House Mouse
FRAGMENT	Thr639-Lys891
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	
MOLECULAR WEIGHT	28.83kDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q7TPD0
ALIAS/SYNONYMS	Int3, SOSS complex subunit A, Sensor of single-strand DNA complex subunit A, SOSS-A, Sensor of ssDNA subunit A