

NACHT, LRR and PYD domains-containing protein 14(Nlrp14)

OVERVIEW

PRODUCT REFERENCE	PX-P4736
ORIGIN SPECIES	
FRAGMENT	Leu401-Phe903
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	68.3kDa
FORM	
BUFFER	0.02% Sarcosyl, PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q6B966
ALIAS/SYNONYMS	Germ cell specific leucine-rich repeat NTPase, NALP-iota,Nalp14